

## Product Datasheet

### MEK3 Antibody (orb1730436)

<b>Description</b>	Dual specificity kinase. Is activated by cytokines and environmental stress in vivo. Catalyzes the c...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, ICC, IHC, IP, WB
<b>Immunogen</b>	A synthesized peptide derived from human MEK3
<b>Preservatives</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 39 kDa; Observed MW: 39 kDa
<b>Uniprot ID</b>	<a href="#">P46734</a>
<b>Entrez</b>	<a href="#">5606</a>
<b>Dilution Range</b>	IP: 1/20; WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/50-1/200; RCM: 1/50-1/100
<b>Expiration Date</b>	12 months from date of receipt.