

Product Datasheet

Histone H2B (yeast) Antibody (orb1730506)

Description	Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessib...
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthesized peptide derived from yeast Histone H2B (yeast)
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Uniprot ID	P02293
Entrez	851810
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.