

Product Datasheet

Phospho-IKK alpha (Thr23) Antibody (orb1730622)

Description	Plays an essential role in the NF-kappa-B signaling pathway which is activated by multiple stimuli s...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide of human CHUK
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Uniprot ID	O15111
Entrez	1147
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.