
Product Datasheet

Phospho-Artemis (Ser516) Antibody (orb1730636)

Description

Required for V(D)J recombination, the

antigen-binding domains of...

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

ELISA, IHC, WB

Immunogen

The antiserum was produced against synthesized peptide derived from human Artemis around the phosphorylation site of Ser516. AA range:482-531

Preservatives

Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

Calculated MW: 78 kDa; Observed MW: 78 kDa

Uniprot ID
[Q96SD1](#)
Entrez
[64421](#)
Dilution Range

WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000

Expiration Date

12 months from date of receipt.