
Product Datasheet

Phospho-FOXO1/3 (Ser322/S325) Antibody (orb1730657)

Description

Transcription factor which acts as a

oxidative stress. In the presenc...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IHC, WB

Immunogen

The antiserum was produced against synthesized peptide derived from human FOXO1A/3A around the phosphorylation site of Ser322 and Ser325. AA range:291-340

Preservatives

Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

Calculated MW: 70 kDa; Observed MW: 97 kDa

Uniprot ID

[O43524](#), [Q12778](#)

Entrez

[2308](#), [2309](#)

Dilution Range

WB: 1/500-1/1000; IHC: 1/50-1/100; ELISA: 1/10000

Expiration Date

12 months from date of receipt.