

Product Datasheet

TriMethyl-Histone H4 (Lys20) Antibody (orb1730703)

Description	The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the D...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human Histone H4 (tri methyl K20)
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	Calculated MW: 11 kDa; Observed MW: 11 kDa
Uniprot ID	P62805
Entrez	121504
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.