

Product Datasheet

Insulin Degrading Enzyme Antibody (orb1730948)

Description	This gene encodes a zinc metallopeptidase that degrades intracellular insulin, and thereby terminate...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	A synthetic peptide of human Insulin degrading enzyme/IDE
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 118 kDa; Observed MW: 118 kDa
Uniprot ID	P14735
Entrez	3416
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.