

Product Datasheet

5T4 Antibody (orb1731129)

Description	5T4 is a 72 kDa transmembrane glycoprotein that exhibits restricted expression in human and mouse ad...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human 5T4
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 46 kDa; Observed MW: 78 kDa
Uniprot ID	Q13641
Entrez	7162
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.