

Product Datasheet

Retinoid X Receptor alpha Antibody (orb1731420)

Description	Receptor for retinoic acid. Retinoic acid receptors bind as heterodimers to their target response el...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	A synthetic peptide of human Retinoid X Receptor alpha
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Uniprot ID	P19793
Entrez	6256
Dilution Range	IP: 1/20; WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.