

Product Datasheet

PRMT1 Antibody (orb1731485)

Description	Arginine methyltransferase that methylates (mono and asymmetric dimethylation) the guanidino nitroge...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human PRMT1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Uniprot ID	Q99873
Entrez	3276
Dilution Range	WB: 1/500-1/1000; IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.