

Product Datasheet

M-CSF Receptor Antibody (orb1731719)

Description	Tyrosine-protein kinase that acts as cell-surface receptor for CSF1 and IL34 and plays an essential ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of mouse MCSF Receptor
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 109 kDa; Observed MW: 109 kDa
Uniprot ID	P09581
Entrez	12978
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.