

Product Datasheet

KAT2A Antibody (orb1731932)

Description	Protein lysine acyltransferase that can act both as a acetyltransferase and succinyltransferase, dep...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant protein of human GCN5
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 94 kDa; Observed MW: 94 kDa
Uniprot ID	Q92830
Entrez	2648
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.