

Product Datasheet

KDM1A Antibody (orb1732592)

Description	Histone demethylase that demethylates both 'Lys-4' (H3K4me) and 'Lys-9' (H3K9me) of histone H3, ther...
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Purified recombinant human KDM1/LSD1 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Calculated MW: 93 kDa; Observed MW: 110 kDa
Uniprot ID	O60341
Entrez	23028
Dilution Range	IP: 1/20; WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.