

Product Datasheet

RBBP5 Antibody (orb1732724)

Description	As part of the MLL1/MLL complex it is involved in methylation and dimethylation at 'Lys-4' of histon...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Purified recombinant human RBBP5 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	Calculated MW: 59 kDa; Observed MW: 70 kDa
Uniprot ID	Q15291
Entrez	5929
Dilution Range	WB: 1/500-1/1000; ICC: 1/50-1/200
Expiration Date	12 months from date of receipt.