

Product Datasheet

PRMT1 Antibody (orb1732794)

Description	Arginine methyltransferase that methylates (mono and asymmetric dimethylation) the guanidino nitroge...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified recombinant human PRMT1 protein fragments expressed in E.coli.
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Uniprot ID	Q99873
Entrez	3276
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.