

Product Datasheet

HDAC4 Antibody (orb1732847)

Description	Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2...
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Purified recombinant human HDAC4 protein fragments expressed in E.coli
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	Calculated MW: 119 kDa; Observed MW: 140 kDa
Uniprot ID	P56524
Entrez	9759
Dilution Range	IP: 1/20; WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.