
Product Datasheet

**Phospho-TAOK1/2/3 (Ser181/Ser181/Ser177)
Antibody (orb1733006)**

Description

Serine/threonine-protein kinase

as p38/MAPK14 stress-activated MA...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IHC, WB

Immunogen

A synthetic phosphopeptide corresponding to residues surrounding Ser181 of human TAOK1

Preservatives

50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin

Concentration

1mg/ml

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

Rabbit IgG

Clonality

Monoclonal

MW

Calculated MW: 116 kDa; Observed MW: 105,116,138 kDa

Uniprot ID
[Q9UL54](#), [Q9H2K8](#), [Q7L7X3](#)
Entrez
[57551](#), [51347](#), [9344](#)
Dilution Range

WB: 1/500-1/1000; IHC: 1/50-1/100

Expiration Date

12 months from date of receipt.