

Product Datasheet

Phospho-CDC6 (Ser106) Antibody (orb1733102)

Description	Involves in the initiation of DNA replication. Also participates in checkpoint controls that ensure ...
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser106 of human Cdc6
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Uniprot ID	Q99741
Entrez	990
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.