

---

## Product Datasheet

### Phospho-Cannabinoid Receptor I (Ser316) Antibody (orb1733106)

**Description**

This gene encodes one of two

cannabinoids, principally delta-9-tetrahydro...

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

WB

**Immunogen**

A synthetic phosphopeptide corresponding to residues surrounding Ser316 of human Cannabinoid Receptor I

**Preservatives**

50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin

**Concentration**

1mg/ml

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Isotype**

Rabbit IgG

**Clonality**

Monoclonal

**MW**

Calculated MW: 53 kDa; Observed MW: 53 kDa

**Uniprot ID**

**P21554**

**Entrez**

**1268**

**Dilution Range**

WB: 1/500-1/1000

**Expiration Date**

12 months from date of receipt.