

Product Datasheet

COL1A1 Antibody (orb1733149)

Description	< span style="font-family:SimHei;">< span style="font-family:"color:#000000;font-size:14px;">This ge...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.03% Proclin 300, pH 7.3.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	139 kDa
Entrez	1277
Dilution Range	IHC: 1/100-1/200
Expiration Date	12 months from date of receipt.