


Product Datasheet

Granzyme B Antibody (orb1733249)

Description	Granzyme B Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Between 197 and C-terminus
Target	Granzyme B
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	250 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P04187
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Detection of mouse Granzyme B by western...