
Product Datasheet

**Mouse FOLR1 Protein, Fc Tag (MALS verified)
(orb1733258)**

Description

Mouse FOLR1 Protein, Fc Tag is

(HEK293). It contains AA Thr 25 - Ser
...**Reactivity**

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

FOLR1

Preservatives

50 mM Tris, 100 mM Glycine, PH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 50.7 kDa. The protein migrates as 60-95 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse FOLR1 Protein, Fc Tag (orb1733258) is expressed from human 293 cells (HEK293). It contains AA Thr 25 - Ser 232 (Accession # P35846-1).

MW

50.7 kDa

Uniprot ID**P35846****Expiration Date**

6 months from date of receipt.