

---

## Product Datasheet

**Mouse PSMA / FOLH1 Protein, Fc Tag (MALS  
verified) (orb1733259)**

**Description**

Mouse PSMA Protein, Fc Tag is

(HEK293). It contains AA Ile 44 - Ala 7...

**Reactivity**

Mouse

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

PSMA

**Preservatives**

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

N-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 106.0 kDa. The protein migrates as > 116 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Mouse PSMA Protein, Fc Tag (orb1733259) is expressed from human 293 cells (HEK293). It contains AA Ile 44 - Ala 752 (Accession # O35409-1).

**MW**

106.0 kDa

**Uniprot ID**[O35409](#)**Expiration Date**

6 months from date of receipt.