

Product Datasheet

Human FGF-4 Protein, His Tag (orb1733261)

Description	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF famil...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	FGF-4
Preservatives	MOPS, Na2SO4, Arginine and EDTA
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 19.8 kDa. The protein migrates as 23-24 kDa under reducing (R) condition (SDS-PAGE).
Purity	90%
Source	Human FGF-4 Protein, His Tag (orb1733261) is expressed from E. coli cells. It contains AA Ser 54 - Leu 206 (Accession # P08620-1).
MW	19.8 kDa
Uniprot ID	P08620
Expiration Date	6 months from date of receipt.