
Product Datasheet

**Mouse TPBG / 5T4 Protein, Fc Tag (MALS verified)
(orb1733266)**

Description

Mouse TPBG Protein, Fc Tag is

(HEK293). It contains AA Ser 32 - Ser 3...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TPBG

Preservatives

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 62.2 kDa. The protein migrates as 80-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse TPBG Protein, Fc Tag (orb1733266) is expressed from human 293 cells (HEK293). It contains AA Ser 32 - Ser 361 (Accession # Q9Z0L0-1).

MW

62.2 kDa

Uniprot ID[Q9Z0L0](https://www.uniprot.org/uniprot/Q9Z0L0)**Expiration Date**

6 months from date of receipt.