
Product Datasheet

**Human PD-1 / PDCD1 Protein, Twin-Strep Tag
(MALS verified) (orb1733267)**

Description

Human PD-1 Protein, Twin-Strep

293 cells (HEK293). It contains AA Leu 25...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

PD-1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-Twin Strep

Note

For research use only

Application notes

This protein carries a twin strep tag at the C-terminus. The protein has a calculated MW of 19.9 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human PD-1 Protein, Twin-Strep Tag (orb1733267) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gln 167 (Accession # Q15116-1).

MW

19.9 kDa

Uniprot ID[Q15116](#)**Expiration Date**

6 months from date of receipt.