
Product Datasheet

**Mouse EGF R Protein, Fc Tag (MALS verified)
(orb1733268)**

Description

Mouse EGF R Protein, Fc Tag is

(HEK293). It contains AA Leu 25 - Ser
...**Reactivity**

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

EGF R

Preservatives50 mM Tris, 100 mM Glycine, 25 mM
Arginine, 150 mM NaCl, pH7.5**Form/Appearance**

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 95.7 kDa. The protein migrates as 110-130 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Mouse EGF R Protein, Fc Tag (orb1733268) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser 647 (Accession # Q01279-1).

MW

95.7 kDa

Uniprot ID[Q01279](#)**Expiration Date**

6 months from date of receipt.