
Product Datasheet

Mouse Mesothelin / MSLN Protein, Fc Tag (orb1733269)

Description

Mouse Mesothelin Protein, Fc Tag is

(HEK293). It contains AA Asp 298 ...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Mesothelin

Preservatives

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 60.6 kDa. The protein migrates as 67-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse Mesothelin Protein, Fc Tag (orb1733269) is expressed from human 293 cells (HEK293). It contains AA Asp 298 - Ser 600 (Accession # Q61468-1).

MW

60.6 kDa

Uniprot ID[Q61468](#)**Expiration Date**

6 months from date of receipt.