
Product Datasheet

**Human MICA (24-202) Protein, Mouse IgG2a Fc Tag
(MALS verified) (orb1733272)**

Description

Human MICA (24-202) Protein,

expressed from human 293 cells (HEK293). It conta...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

MICA

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-mFc(mlgG2a)

Note

For research use only

Application notes

This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 47.8 kDa. The protein migrates as 57-64 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human MICA (24-202) Protein, Mouse IgG2a Fc Tag (orb1733272) is expressed from human 293 cells (HEK293). It contains AA Glu 24 - Arg 202 (Accession # Q29983-1).

MW

47.8 kDa

Uniprot ID**Q29983****Expiration Date**

6 months from date of receipt.