

Product Datasheet

Mouse Beta-NGF Protein, Tag Free (MALS verified) (orb1733273)

Description	Mouse Beta-NGF Protein, Tag Free is expressed from human 293 cells (HEK293). It contains AA Ser 122 ...
Reactivity	Mouse
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	Beta-NGF
Preservatives	20 mM NaAC, 150 mM NaCl, pH5.5
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 13.5 kDa. The protein migrates as 14 kDa under reducing (R) condition (SDS-PAGE).
Purity	95%
Source	Mouse Beta-NGF Protein, Tag Free (orb1733273) is expressed from human 293 cells (HEK293). It contains AA Ser 122 - Gly 241 (Accession # P01139).
MW	13.5 kDa
Uniprot ID	P01139
Expiration Date	6 months from date of receipt.