
Product Datasheet

**Rat TPBG / 5T4 Protein, His Tag (MALS verified)
(orb1733279)**

Description

Rat TPBG Protein, His Tag is

(HEK293). It contains AA Ser 32 - Ser 36...

Reactivity	Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TPBG
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.7 kDa. The protein migrates as 50-70 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Rat TPBG Protein, His Tag (orb1733279) is expressed from human 293 cells (HEK293). It contains AA Ser 32 - Ser 361 (Accession # Q5PQV5).
MW	37.7 kDa
Uniprot ID	Q5PQV5
Expiration Date	6 months from date of receipt.