
Product Datasheet

Human LRP-1 / cluster IV Protein, Fc Tag (MALS verified) (orb1733282)

Description

Human cluster IV Protein, Fc Tag is

(HEK293). It contains AA Ser 333...

Reactivity

Human

Endotoxins

0.1 EU per µg

Conjugation

Unconjugated

Target

LRP-1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 76.6 kDa. The protein migrates as 80-95 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human cluster IV Protein, Fc Tag (orb1733282) is expressed from human 293 cells (HEK293). It contains AA Ser 3332 - Asp 3779 (Accession # Q07954-1).

MW

76.6 kDa

Uniprot ID

[Q07954](#)

Expiration Date

6 months from date of receipt.