
Product Datasheet

Human LRP-1 / cluster III Protein, Fc Tag (MALS verified) (orb1733283)

Description

Human cluster III Protein, Fc Tag is

(HEK293). It contains AA Ser 252...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LRP-1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 72.6 kDa. The protein migrates as 80-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human cluster III Protein, Fc Tag (orb1733283) is expressed from human 293 cells (HEK293). It contains AA Ser 2522 - Ile 2941 (Accession # Q07954-1).

MW

72.6 kDa

Uniprot ID[Q07954](#)**Expiration Date**

6 months from date of receipt.