

Product Datasheet

Anti-dsRNA Comparison Set (orb1733300)

Description	SCICONS K1 monoclonal antibody recognises dsRNA with similar affinity to the widely used SCICONS J2 ...
Conjugation	Unconjugated
Form/Appearance	lyophilised
Storage	store at -20 °C to -80 °C upon reconstitution for long-term storage. avoid freeze/thaw cycles, store in aliquot. After adding 10 mM sodium azide the undiluted antibody (1 µg/µl) can be stored at 4°C for a short period of time
Note	For research use only
Solubility (25°C)	To prepare a 1 µg/µl PBS antibody solution add 100 µl, sterile DNase/RNase-free distilled water. As a result of the lyophilisation procedure, the reconstituted antibody may contain small amounts of denatured protein in the form of aggregates that may interfere with some applications such as immunohistochemistry, (e.g. by giving high backgrounds). We therefore highly recommend to spin down (microcentrifuge) the reconstituted antibody before use, and to use the supernatant only.
Expiration Date	12 months from date of receipt.