

Product Datasheet

HighFidelity IR750 PCR Labeling Kit (orb1733324)

| | |
|--------------------------|--|
| Description | HighFidelity IR750 PCR Labeling Kit is designed to produce randomly IR750-modified DNA probes by PCR... |
| Storage | store at -20 °C. avoid freeze/thaw cycles, store dark |
| Note | For research use only |
| Application notes | Spectroscopic Propertie: λ_{exc} 753 nm, λ_{em} 775 nm, ϵ 255.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5). |
| Expiration Date | 12 months from date of receipt. |