

## Product Datasheet

### Furin (orb1733332)

<b>Description</b>	Furin (proprotein convertase subtilisin/kexin type 3, PCSK3, PACE) belongs to the proprotein convert...
<b>Form/Appearance</b>	liquid (supplied in 10 mM Hepes/NaOH pH 7.4, 100 mM NaCl, 2mM CaCl <sub>2</sub> )
<b>Storage</b>	store at -80 °C. avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Purity</b>	> 90% (SDS-PAGE)
<b>Source</b>	human cell line HEK293S
<b>MW</b>	Theoretical MW: 52.3 kDa
<b>Expiration Date</b>	6 months from date of receipt.