

Product Datasheet

Legumain (orb1733562)

Description	Legumain is a cysteine protease with a strict specificity for cleaving after asparagine and, to a le...
Form/Appearance	liquid (supplied in 20 mM citric acid pH 4.0, 100 mM NaCl, 2 mM DTT)
Concentration	500 µg/ml
Storage	store at -20 °C. avoid freeze/thaw cycles
Note	For research use only
Clonality	Recombinant
Purity	> 90% (SDS-PAGE)
Source	Leishmania
MW	Theoretical MW: 36 kDa
Expiration Date	6 months from date of receipt.