

## Product Datasheet

### SAP-Exo Kit (orb1733626)

<b>Description</b>	The Kit contains two hydrolytic enzymes, recombinant Shrimp Alka-line Phosphatase (rSAP) and Exonucl...
<b>Storage</b>	store at -20 °C. avoid freeze/thaw cycles
<b>Note</b>	For research use only
<b>Application notes</b>	SAP-Exo Kit is a quick, easy and reliable enzymatic cleanup reagent for PCR product cleanup. Removes excess primers and dNTPs; Fast 20 minute protocol; Add directly to PCR product; 100% Sample Recovery; Scalable for different reaction sizes; No interference on downstream applications; Easy to automate.
<b>Expiration Date</b>	12 months from date of receipt.