

Product Datasheet

Total RNA Purification Kit - Column Kit (orb1733629)

| | |
|------------------------|---|
| Description | Total RNA Purification Kit is designed for rapid, high purity and high yield isolation of total RNA from small amounts of various samples including blood, animal and plant tissue, bacteria and viruses. The spin column based method allows complete removal of inhibitors such as divalent cations and proteins. Due to elimination of phenol, handling of the kit is safe and no harmful waste is produced. The purified total RNA can be used in a number of downstream applications. The kit allows the purification of up to 100 µg RNA per preparation. |
| Storage | store at ambient temperature |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)