

Product Datasheet

Anchored Oligo (dT)20 (orb1733638)

Description	Anchored Oligo (dT) Primer hybridizes with highest efficiency to the poly(A) tail of mRNA and is recommended as primer for first stand cDNA synthesis, cDNA labeling or RT-PCR applications.
Form/Appearance	clear aqueous solution
Concentration	100 μ M
Storage	store at -20 °C. avoid freeze/thaw cycles
Note	For research use only
Purity	by HPLC
MW	Theoretical MW: 6639 g/mol
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)