

## Product Datasheet

# SCRIPT Direct RT-qPCR ProbesMaster highROX (orb1733650)

|                        |   |
|------------------------|---|
| <b>Description</b>     | SCRIPT Direct RT-qPCR ProbesMaster highROX is designed for quantitative real-time analysis of target... |
| <b>Form/Appearance</b> | liquid  |
| <b>Concentration</b>   | 2x  |
| <b>Storage</b>         | store at -20 °C. avoid freeze/thaw cycles. stable at 4 °C for up to 4 weeks                             |
| <b>Note</b>            | For research use only   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |