

Product Datasheet

SCRIPT RT-qPCR SybrMaster highROX (orb1733657)

Description	SCRIPT RT-qPCR SybrMaster highROX is designed for quantitative real-time analyses of RNA templates u...
Form/Appearance	liquid
Concentration	2x conc
Storage	store at -20 °C. avoid freeze/thaw cycles, store dark. stable at 4 °C for up to 4 weeks
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 494 nm (bound to DNA), λ_{em} 521 nm, (bound to DNA).
Expiration Date	12 months from date of receipt.