

## Product Datasheet

### Direct qPCR ProbesMaster (orb1733665)

**Description**

Direct qPCR ProbesMaster is designed for quantitative real-time analysis of target DNA directly from blood, swabs and animal-or plant tissue. The mix allows robust amplification avoiding the requirement of any prior DNA purification procedures. The mix is recommended for use with dual-labeled fluorescent probes, e.g. TaqMan, Molecular Beacon or Scorpion probes. It provides a powerful tool for multiplex-quantification of sample DNA in a broad dynamic range with exceptional sensitivity and precision. The mix contains all reagents required for qPCR (except template, primer and labeled fluorescent probe) in a premixed 2 x concentrated ready-to-use solution. High robustness, reliability and sensitivity of the mix are based on a an antibody-blocked hot start polymerase in combination with an optimized and well-balanced buffer system. The mix ensures fast and easy preparation with a minimum of pipetting steps and is specially recommended for: Direct detection of viral or bacterial DNA in nasal or throat swabs; Direct PCR from blood samples; Direct amplification of target DNA from various tissue samples; Point-of-Care diagnostics.

**Form/Appearance**

liquid

**Concentration**

2x

**Storage**

store at -20 °C. avoid freeze/thaw cycles,. stable at 4 °C for up to 4 weeks

**Note**

For research use only

**Expiration Date**

12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)