

## Product Datasheet

# Saphir Bst Turbo GreenMaster Lyophilisate (orb1733666)

<b>Description</b>	Saphir Bst Turbo GreenMaster Lyophilisate is designed for isothermal amplification of DNA. The mix i...
<b>Storage</b>	store at ambient temperature. Store in an aluminium-coated bag or on a dry place. Lyophilisates may hydrate at humidity levels > 70 % when sealing is opened
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.