

## Product Datasheet

### qPCR ProbesMaster lowROX (orb1733687)

<b>Description</b>	qPCR ProbesMaster lowROX is designed for quantitative real-time analysis of DNA samples using DNA pr...
<b>Form/Appearance</b>	liquid
<b>Concentration</b>	2x conc
<b>Storage</b>	store at -20 °C. avoid freeze/thaw cycles, store dark. Storage at 4 °C for up to 3 months possible
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.