

Product Datasheet

PCR-grade Water (orb1733694)

Description	PCR-grade Water is intended for use in molecular biology applications including PCR and RT-PCR. The ultra-pure and sterile filtered water is manufactured free of detectable inhibitors, contaminants or enzymatic activity.
Form/Appearance	liquid
Storage	store at 4 °C or -20 °C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)