

## Product Datasheet

# SCRIPT RT-qPCR ProbesMaster Lyophilisate (orb1733721)

<b>Description</b>	SCRIPT RT-qPCR ProbesMaster Lyophilisate is designed for performing maximum sensitive and specific R...
<b>Storage</b>	store at ambient temperature. Store in an aluminium-coated bag or on a dry place. Lyophilisates may hydrate at humidity levels > 70 % when sealing is opened
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.