

Product Datasheet

Fluorescent 500 bp DNA Ladder (orb1734244)

Description	Fluorescent 500 bp DNA Ladder consist of 10 DNA fragments ranging from 500 bp to 5 kb. The prestaine...
Form/Appearance	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Green, Bromophenol Blue and Xylene Cyanol FF; Color: blue
Concentration	105 ng/μl
Storage	store at -20 °C. store dark. Short term storage (up to 1 month) at 4 °C possible
Note	For research use only
Application notes	Recommended Load: 5 μl per lane.
Expiration Date	12 months from date of receipt.