

## Product Datasheet

### Fluorescent 200 bp DNA Ladder (orb1734245)

<b>Description</b>	Fluorescent 200 bp DNA Ladder consist of 10 DNA fragments ranging from 200 bp to 2 kb. The prestaine...
<b>Form/Appearance</b>	liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Green, Bromophenol Blue and Xylene Cyanol FF; Color: blue
<b>Concentration</b>	105 ng/μl
<b>Storage</b>	store at -20 °C. store dark. Short term storage (up to 1 month) at 4 °C possible
<b>Note</b>	For research use only
<b>Application notes</b>	<b>Recommended Load:</b> 5 μl per lane.
<b>Expiration Date</b>	12 months from date of receipt.