

Product Datasheet

Fluorescent 100 bp DNA Ladder (orb1734246)

| | |
|--------------------------|--|
| Description | Fluorescent 100 bp DNA Ladder consist of 10 DNA fragments ranging from 100 bp to 1 kb. The prestaine... |
| Form/Appearance | liquid, DNA fragment mix in 10 mM Tris-HCl pH 7.5, 10 mM EDTA, 10% (w/v) Glycerol, SYBR® Tartrazine and Xylene Cyanol FF; Color: green |
| Concentration | 105 ng/μl |
| Storage | store at -20 °C. store dark. Short term storage (up to 1 month) at 4 °C possible |
| Note | For research use only |
| Application notes | Recommended Load: 5 μl per lane. |
| Expiration Date | 12 months from date of receipt. |